


Certificate

Aidian Oy

using

QuikRead go HbA1c

participated in the IFCC HbA1c Certification Programme to demonstrate traceability to the IFCC Reference Measurement Procedure and performed as shown below.




Total Error =	1.7	mmol/mol
Bias =	-0.1	mmol/mol
Imprecision =	1.6	%
Grade =	Silver	

Criteria derived from the IFCC model for Quality Targets HbA1c (Clin Chem 2015;61 : 752-59)

Date of Certification : 01 June 2023

Date of Expiry : 01 June 2024



IFCC Network Coordinator
C. Siebelder